


NASAL ASPIRATORS

Nasal Aspirators are most commonly used to remove mucus from babies' noses but can be used in a variety of other applications as well. The tips on this product can be easily removed for cleaning. Our bulbs are made from medical-grade Sarolit™ and are not made with BPA, phthalates or latex. Our products are also made in one piece using the rotational molding process so that there are no bonded seams, which could fail during service or foster bacterial buildup.

This product is currently offered in the standard colors listed in the table. It can be made in any of the colors listed below or a custom color at customer request. Products can also be printed with customer logos or otherwise modified to meet customer specifications.

- 01- white
- 02- pink
- 03- light blue
- 04- black
- 05- red
- 06- kelly green
- 07- purple
- 08- cobalt blue
- 12- oyster gray
- 13- medium gray



Visit the Albert Manufacturing USA website for all accessory specs.

Small Aspirator with Curved Tip

Name and Description	Volume	Length (L)	Diameter (D)	Tip Material	Part Number*	Standard Bulb Colors
1oz Nasal Aspirator with Curved Soft Tip, Sarolit™	1oz	87mm	41mm	Sarolit™	CN/N01Mxx-20N-20	01 03



Aspirator with Tapered Tip

Name and Description	Volume	Length (L)	Diameter (D)	Tip Material	Part Number*	Standard Bulb Colors
2oz Nasal Aspirator with Tapered Hard Tip, PC	2oz	88mm	41mm	Polycarbonate	DE/N02Mxx-24W-21	01 03
2oz Nasal Aspirator with Tapered Hard Tip, PP	2oz	88mm	41mm	Polypropylene	DE/N02Mxx-24W-38	01 03
2oz Nasal Aspirator with Tapered Soft Tip, Sarolit™	2oz	88mm	41mm	Sarolit™	DE/N02Mxx-24W-37	01 03



Aspirator with Standard Tip

Name and Description	Volume	Length (L)	Diameter (D)	Tip Material	Part Number*	Standard Bulb Colors
2oz Nasal Aspirator with Standard Hard Tip, PC	2oz	88mm	41mm	Polycarbonate	DE/N02Mxx-24W-36	01 03
2oz Nasal Aspirator with Standard Hard Tip, PP	2oz	88mm	41mm	Polypropylene	DE/N02Mxx-24W-25	01 03
2oz Nasal Aspirator with Standard Soft Tip, Sarolit™	2oz	88mm	41mm	Sarolit™	DE/N02Mxx-24W-35	01 03



* The "xx" in the Part Number is a placeholder for the color code as shown above