

# STRESS RELIEF BALLS

The stress relief balls are made with soft and pleasant bumps across the surface for body stimulation. They can be used to ease tension in muscles or to aid in rehabilitation and mobilization of patients. The stress relief balls are made from medical-grade Sarolite™ and are not made with BPA, phthalates or latex. Our products are also made in one piece using the rotational molding process so that there are no bonded seams, which could fail during service or foster bacterial buildup.

This product is currently offered in the standard colors listed in the table. It can be made in any of the colors listed below or a custom color at customer request. Products can also be printed with customer logos or otherwise modified to meet customer specifications.

This product is made in our German plant.

- 01- white
- 02- pink
- 03- light blue
- 04- black
- 05- red
- 06- kelly green
- 07- purple
- 08- cobalt blue
- 12- oyster gray
- 13- medium gray

Visit the Albert Manufacturing USA website for all accessory specs.

Name and Description	Size	Part Number*	Standard Colors
5.5cm Stress Relief Ball	5.5cm	9005.2.xx000	05 06 08
6.0cm Stress Relief Ball	6.0cm	9006.2.xx000	05 06 08
7.0cm Stress Relief Ball	7.0cm	9007.2.xx000	05 06 08

